

# danumed® Extension Set for Gastrostomy Button


*Extension set for delivery of enteral nutrition, liquid medication and decompression.*

*Innovative features optimise the flow rate and simplify the handling for the user.*

*Cost-efficient solution with usage for up to 14 days.*



## COMFORT



### Rotating Connector




- Allows for the easy connection and release of ENFit® connections.
- Even very tight connections can be easily released.
- Only the rotating Swivel is turned and not the entire connector. This feature **avoids any twisting of the enteral tube** which is both painful for the patient and inconvenient for the user.
- All danumed ENFit® connectors offer the advantages of the ENSwivel®: **ENFit®**, **Y-ENFit®** and **Bolus-ENFit®**.




### High-End Button Connector

- The locking mechanism with optimised design ensures an **easy and tight connection** to the button.



The rounded kink in the inner lumen optimises the flow rate, independent of the tube lumen.



Ergonomic form of the button connector with ribbed profile for a better grip.

## SAFETY

### Increased safety for the patient

- The ENFit® connector according to DIN EN ISO 80369-3 **avoids tubing misconnections** and increases patient safety.
- All ENFit® connectors from danumed are made of **TRITAN™**, a plastic material resistant to aggressive solutions/medicines.
- All extension sets with ENFit® connector are **free from latex, PVC, DEHP and BPA**.
- The biocompatible tube made of polyurethane is **kink resistant** and **reduces the hazard of tube occlusions** and the number of alarm signals of the enteral pump.



## danumed® Extension Set with angled Button Connector



REF	Description	Length (cm)	Feeding connector	Latex-free	PVC-free	DEHP-free	BPA-free
DA 1011 3011	Extension Set <b>sterile</b> , 30 cm, ENFit®, ENSwivel®	30 cm	ENFit®	•	•	•	•
DA 1011 6011	Extension Set <b>sterile</b> , 60 cm, ENFit®, ENSwivel®	60 cm	ENFit®	•	•	•	•
DA 1011 0521	Extension Set <b>sterile</b> , 5 cm, Y-ENFit®, ENSwivel®	5 cm	Y-ENFit®	•	•	•	•
DA 1011 3021	Extension Set <b>sterile</b> , 30 cm, Y-ENFit®, ENSwivel®	30 cm	Y-ENFit®	•	•	•	•
DA 1011 6021	Extension Set <b>sterile</b> , 60 cm, Y-ENFit®, ENSwivel®	60 cm	Y-ENFit®	•	•	•	•

Delivered with right-angled button connector and **standard bore tube**. Individually packaged. 10 pieces in a box, 10 boxes in a shipping box.



## danumed® Extension Set with straight Button Connector



REF	Description	Length (cm)	Feeding connector	Latex-free	PVC-free	DEHP-free	BPA-free
DA 1011 3032	Extension Set <b>sterile</b> , 30 cm, Bolus-ENFit®, ENSwivel®	30 cm	Bolus-ENFit®	•	•	•	•
DA 1011 6032	Extension Set <b>sterile</b> , 60 cm, Bolus-ENFit®, ENSwivel®	60 cm	Bolus-ENFit®	•	•	•	•

Delivered with straight button connector and **wide bore tube**. Individually packaged. 10 pieces in a box, 10 boxes in a shipping box. The Bolus Extension Sets are available on request with funnel connector instead of ENFit® connector.

### Related Products to danumed® Extension Sets:



danuButton®  
Gastrostomy Button



INNOVATORS IN ENTERAL CARE

For more information, please contact us on (02) 9882 3666 or at [amsl@amsl.com.au](mailto:amsl@amsl.com.au)

[amsl.com.au](http://amsl.com.au)

AMSL is a subsidiary of Dexcom. PR-100-582 June 2022

